


BBMPulser 4B

user manual



Welcome

BBMPulser 4B generates pulsed magnetic fields which create micro-currents of electricity that work with the body's natural electricity for general health and well-being.

The magnetic pulser is designed to be used independently, but can also be used together with Silver Pulser, and Ozone Generator as per dr. Robert Becks Protocol titled "Take back your power"

Before Using

Please read this manual thoroughly before using your magnetic pulser. For the most up-to-date manual, please visit our website [**www.bbmpulser.org/manuals**](http://www.bbmpulser.org/manuals)

We have made every effort to ensure the information in your manual is up-to-date at the time of printing ... since that time however, new information may have become available.

Learn more

To learn more about the Magnetic Pulsers we invite you to visit the BBMPulser website: [**www.bbmpulser.org**](http://www.bbmpulser.org) The website offers Gallery of pictures, Blog, as well as Frequently Asked Questions (coming soon), The Magnetic Pulsers Products User Guide, and more.

Complete Unit includes

The unit Magnetic Pulser model BBMPulser 4B

120VAC Power cord

Product user manual

Carry bag

- 24 gram 128 aluminum foil sheets for bend and repel tests

Do not Use

1. Do not use with a pacemaker.
2. Do not use in the following situations as safety has not been explored:
 - a. in the case of pregnancy.
 - b. in the area of an implanted device such as a shunt, etc.
 - c. in the area of external or internal bleeding wounds, including menstruation.
 - d. in the area of a metal implant. The metals may get pull or repel by magnetic force
3. As with all electrical products do not use near water in high humidity or high heat rooms (sauna), when driving a car or when operating heavy equipment. Try NOT to fall asleep with the unit turned on, the coil will get hot and may cause burns if left in one spot.

Cautions

1. If you experience nerve pain or an increase in nerve pain while using the unit, please discontinue using in the area of pain.
2. Caution is advised when using heat-producing or heat-sensitive ointments, lotions or medications as a skin burn may occur.
3. With normal use, the hand paddle will become very warm to hot. There is a 50C thermal switch build in near the coil, when it senses appropriate heat it will place the unit in the stand by mode, red "Protec" LED light upfront will start blinking. Please allow 20 min for the coil to cool down on its own after

which, the Unit will turn On. (Safety future)

4. To prevent damage or loss of information, keep away from homeopathic remedies, credit cards, security cards, audio and video tapes. Also from electronic equipment— computer disks, hearing aids, cell phones, answering machines, Car keys and others valuable electronics.

5. While this technology is generally considered safe, there exists potential for rare individual reactions that cannot reasonably be foreseen. Therefore, your use of the BBMPulser products constitutes your agreement that you are responsible for your decision to use the technology. BBMPulser assumes no responsibility for any direct or indirect loss or damage believed to arise from using the BBMPulser 4B Unit.

6. Do not cover the unit, it requires continues un-restricted airflow for normal operation. In case airflow gets restricted - build in 50C thermal switch will trigger, placing the unit in the “stand by” mode until the unit cools down, disabling power to the coil. Ones unit cools off the thermal switch will automatically reset turning the unit back on. (Safety future)

7. When testing magnetic output with the aluminum foil 24 gram test peace DO NOT place your face or any parts of your body in front of it, when it gets repel with magnetic force it can fly up in the face and cause damage. Perform repel test only outside with unlimited ceiling space.

Summary of Lights

Green Power LED: indicates the unit is on.

Red Blinking: indicates the Unit is in the “stand by” mode and cooling off, it happened because either coil or unit is overheating. Allow 15-20min to cool off and try again.

Blue Blinking: indicates pulse is being released.

Red South Pole indicator "S" molded in the head of the peddle indicates which side of the coil is South, use North side of the coil in 95% of the cases. Meaning the South Pole LED Indicator should be pointing away from the body during use.

How to Use the Magnetic Pulser

1. Insert power plug into an electrical outlet 120VAC
2. Push the ON/OFF button upfront

You will hear the click ones every 4-5 seconds as the Magnetic pulse released, it is normal.

4. The head of the hand paddle should be located on or as close to the body as possible. It does not need to touch the skin as clothing will not interfere with the magnetic field. In some cases for deeper penetration the head of the paddle could be slightly pushed against the body for deeper penetration. The strongest magnetic field is at the surface of the coil.

Note: The pulsed magnetic field from this unit will penetrate more than 20 inches into the body. The bacteria, virus killing zone is around 2-4 inches off the surface of the coil

5. To turn the unit off, push the "Power On/Off" button at any time, you don't need to have unit plugged in while cooling off.

How it should feel:

Most individuals will not feel any sensation while using the magnetic pulser some individuals may feel small to medium electric shock in some places (side of the neck and under arms) most of the body you won't feel anything.

With normal use, the hand paddle will become very warm to hot. When it happens build in thermal switch will place the unit in stand by mode to cool

the coil off

Typical use

Length of Session when used for tissue recovery You can use this Unit for intense deep stimulations or to treat infections. If the purpose is to stimulate then start with couple min and build your way up to an hour/day adding couple min /session/each day

Length of Session when used to deal with infections. Depending on the size and origin of the infection the length of Session can vary between 1 min to an hour/ day.

Length of Use: In some cases all you need is couple of pulses to sterilize the cut or to disinfect something, in other cases you may end up using it 5-6 times/day 20min each time.

For specific purposes, longer and more frequent sessions may be more beneficial. Generally, the Magnetic Pulser can be used as needed.

Placement:

You get cone shape killing zone coming out from both sides of the coil with each pulse released. The cone is 3" wide and 4" tall at the tip.

You only need a few Pulses/area, the unit is putting out strong pulses capable of killing virus or bacteria instantly.

Drink plenty of water: It is important to drink enough water when using the unit. Drinking Distilled and Ozonated water is especially helpful.

Stop Gradually: Reduce the amount of time each day for the last week of use.

Electrical specifications

Your Magnetic Pulser outputs a time-variant pulsed DC magnetic field, 3 milliseconds in length.

Power Requirements:

120Volts AC, peaks at 300Watts **Caution:** Exceeding the voltage rating may result in blown fuse. Please only use 120VAC

Output specifications

Magnetic Field Intensity:

14,000 Gauss (1.4 Tesla) at face of the hand paddle

Pulse Duration: 3 mS (milliseconds).

Weight of the copper in the coil - 400 Gram

Coil Wire diameter: 1.7mm (14AWG)

Unit Peddle made using: Urethane Casting Resin and Flex-it 90

Capacitor Bank: 2,300uf

Discharge Voltages: 550V

Testing the strength of Magnetic Pulse

24 Gram Aluminum foil test piece is included for testing purposes. North Side of the peddle is used for the test. The coil inside the peddle intentionally placed closer to the north facing side. Reason? Because 95% of the time North side of the coil is used, so the coil positioned closer to the North side of the peddle to deliver stronger pulse from the North side.

You can test the force of magnetic pulse from the head of the paddle by placing stack of aluminum foils (included) on the top of the coil and trying on different power settings. BE CAREFUL as test piece will repel in the air 10'

or more. Perform the test outdoors only. Make sure no cars near by, as it may get blown by the air and land on the car causing damage. Do not put your head or any parts of your body in direct line where foil would fly to avoid getting injured. Its recommended to wear safety glasses performing this test.

Troubleshooting:

If you experience issues with the unit please reach us at support@bbmpulser.org

Warranty

The magnetic pulser model 4B is covered by a two (2) year limited warranty. Please keep proof of purchase, Warranty is void if the failure is due to abuse or negligence.

Please Note

Some older houses do not have ground connector. In this case use and adapter allowing to bypass ground terminal. Unit will operate normally with out being grounded Please Use ONLY 120VAC as a power source. Internal circuitry is not designed to work with higher or lower voltages.

Please memorize to use North Side of the coil in 95% of the cases, (South side of the coil is used when trying to grow something on your body for example loss hair, or nails or in some cases the tooth, for any other cases you would use NORTH side of the coil).

Repairs

Please don't attempt to repair the unit, DANGER OF HIGH VOLTAGE exist. Opening up the unit will void your warranty.

If the unit malfunctions, please contact our support for RMA instructions by sending eMail to: support@bbmpulser.org or call 1.425.312.3067 USA.