

USER MANUAL

5 Tesla BBMPulser 6B

Welcome

BBMPulser generates Pulsed Magnetic Fields which create micro-currents of electricity that work with the body's natural electricity for general health and well-being. When used on higher settings (75%, 100%) the BBMPulser 6B generates stronger currents capable of killing Bacteria and Viruses near the surface of the coil.

The Magnetic Pulser is designed to be used independently, but can also be used together with Silver Pulser, and Ozone Generator as per dr. Robert Becks Protocol titled "THE BECK PROTOCOL A FIRST AID KIT OF THE FUTURE"

<https://bbmpulser.org/2019/12/23/the-beck-protocol-a-first-aid-kit-of-the-future/>

Before Using

Please read this manual thoroughly before using your Magnetic Pulser. For the most up-to-date manual, please visit our website www.bbmpulser.org/manuals

Complete Unit includes:

- The unit Magnetic Pulser model BBMPulser 6B 5 Tesla (50,000 Gauss)
- Product user manual

Do not Use

1. Do not use with a pacemaker.
2. Do not use in the following situations as safety has not been explored:
 - a. in the case of pregnancy.
 - b. in the area of an implanted device such as a shunt, etc.
 - c. in the area of external or internal bleeding wounds, including menstruation.
 - d. in the area of a metal implant. The metals containing iron will be pulled to the Coil (incl stainless steel), all other metals will be repelled from the Coil North and South sides by strong Magnetic Fields generated by the Coil.
3. As with all electrical products do not use near water in high humidity or high heat rooms (sauna).

Try NOT to fall asleep with the unit turned on, the Coil will get hot and may cause burns if left in one spot.

Cautions:

1. If you experience nerve pain or an increase in nerve pain while using the unit, please discontinue using in the area of pain.
2. Caution is advised when using heat-producing or heat-sensitive ointments, lotions or medications as a skin burn may occur.
3. With normal use, the hand paddle will become very warm to hot. There is a 60C thermal switch build in near the Coil, when it senses appropriate heat it will trigger Sound Alarm and Red Protec LED, when it happens, please turn off your Unit and allow 30-45min for Coil to cool off. Do not operate unit while sound Alarm is ON.
4. To prevent damage or loss of information, keep away from homeopathic remedies, credit cards, security cards, audio and video tapes. Also from electronic equipment—computer disks, hearing aids, cell phones, answering machines, Car keys and others valuable electronics.
5. While this technology is generally considered safe, there exists potential for rare individual reactions that cannot reasonably be foreseen. Therefore, your use of the BBMPulser products constitutes your agreement that you are responsible for your decision to use the technology. BBMPulser assumes no responsibility for any direct or indirect loss or damage believed to arise from using the BBMPulser 6B Unit.
6. Do not cover the unit, it requires continues unrestricted airflow for normal operation.

In case airflow gets restricted - build in thermal switch will trigger, triggering "sound Alarm", if this happens, please turn off the Unit. Once unit cools off the thermal switch will automatically reset turning the unit back on. (Safety feature)

Summary of Lights:

1. Green Power LED - indicates the Unit is on.
2. Red Blinking LED - indicates the Unit or the Coil is in the "Thermal Protection Alarm" mode and need to cool off, if it happens allow 30-45 min to cool off and try again.
- Blue Blinking LED - indicating pulse is being released.

Your Coil has North and South Sides. North side is the hard side of the Coil (white or red)

In 95% of the cases, North side is used, North side to the body

How to Use the Magnetic Pulser

1. Insert power plug into an electrical outlet 120VAC
2. Make sure your Coil is not sits on any metal, place it on the none metallic surface
3. Switch the ON/OFF switch.
4. Make your power selection, you can select from 2 Tesla to 5 Tesla intensity with rotary switch in front of the unit.

The click will be heard as the pulse is released especially on higher settings, it is normal.

5. The head of the hand paddle should be located on or as close to the body as possible.

It does not need to touch the skin as clothing will not interfere with the magnetic field.

In some cases for deeper penetration the head of the paddle could be slightly pushed against the body during treatment. The strongest magnetic field is at the surface of the Coil

Note: The pulsed magnetic field from this unit will penetrate more then 20 inches into the body.

5. To turn the unit off, Switch the "Power On/Off" at any time, you don't need to unplugged Unit from 120VAC while cooling off.

How it should feel:

When used on high setting you may feel small to medium electric shock in some places, in some areas your muscles will twitch with each pulse

With normal use, the hand paddle will become very warm to hot. When it happens build in thermal switch will trigger, please turn off the unit and allow Coil to cool off.

Typical use length of Session is around 7 min. It is best to begin slowly (3-5 minutes) and then gradually increase the length of each session. Daily use varies from 1 minute (in some cases) up to two hours or more.

Lower settings (2 Tesla) good for deep stimulation, injury recovery, top skin infections

Stronger Settings (3-5 Tesla) good when dealing with infections an when there is a need to reach deeper, the higher on setting you go the deeper it can penetrate.

If the purpose is to stimulate then start with couple min and build your way up to an hour/day adding couple min each session each day

Length of Session when used at full power. This setting is used when dealing with infections.

Depending on the size and origin of the infection the length of session can vary between 1 min to an hour.

Length of Use: In some cases all you need is couple min session to sterilize the cut or to disinfect something, in other cases you may end up using it 5-6 times/day 7 min each time.

For specific purposes, longer and more frequent sessions may be more beneficial.

Generally, the Magnetic Pulser can be used as needed.

Placement: (50% power) The hand paddle may be kept on one location for up to 5 min . It may also be moved at any time to a different location

At 75% and 100% power the hand paddle should not be kept at one location, you should move it around as it pulses.

You will get cone shape of the magnetic flux coming out from both sides of the Coil.

Drink plenty of water: It is important to drink enough water when using the unit.

Drinking Purified and Ozonated water is especially helpful.

Stop Gradually: Reduce the amount of time each day for the last week of use.

Electrical specifications:

Your Magnetic Pulser outputs a time-variant pulsed DC magnetic field, 1-2 milliseconds in length.

Power Requirements 120Volts AC it will peak at 25Amps for short period of time after each pulse and on average it will consume around 400Watts

Caution: Exceeding the voltage rating will result in build in Circuit breaker triggering or it may result in damage to the Unit. Please Do NOT exceed rated Voltage

Output specifications

Magnetic Field Intensity:

@ 50% ~20,000 Gauss at face of the hand paddle. (2 Tesla)

@ 75% ~35,000 Gauss at face of the hand paddle. (3.5 Tesla)

@ 100% 50,000 Gauss (5 Tesla) at face of the hand paddle

Pulse Duration: 1-2 ms (milliseconds).

Weight of the copper in the coil - 500 Gram

Coil wire specs: Copper Strip 15mm wide x 0.5mm thick

Coil ID 1", OD 3"

Unit Peddle made using: 2 components, hard and soft resins Both components are made of Urethane Casting Resin "Flex-it 90" and white or black hard plastic.

Capacitor Bank: 1,350uF Max charge voltage 1.5kV

Unit weight: ~ 70Lb

Warranty:

The Magnetic Pulser model 6 is covered by two (2) year limited warranty. Please keep proof of purchase, Warranty is void if the failure is due to abuse or negligence.

Please Note: Some older houses do not have ground connector. In this case use an adapter allowing to bypass ground terminal. Unit will operate normally without being grounded.

Please Use ONLY 120VAC as a power source. Internal circuitry is not designed to work with higher or lower voltages.

Please memorize to use North Side of the coil in 95% of the cases, (South side of the coil is used when trying to grow something on your body for example loss hair, or nails or in some cases the tooth, for any other reasons you would use NORTH side of the coil).

Repairs:

Please don't attempt to repair the unit, DANGER OF HIGH VOLTAGE exist. Opening up the unit will void your warranty and will put you in high risk of getting shocked by High Voltage.

If the unit malfunctions, please contact our support for RMA instructions by sending eMail to:

You can reach us at: support@bbmpulser.org USA.